

CAB2302

Product Information

Product SKU: CAB2302	Gene ID: 7827	Size: 20uL, 100uL
Clone No: ARC1907	Host Species: Rabbit	Reactivity: Human

Additional Information

Observed MW: 35kDa	Conjugate: Unconjugated
Calculated MW: 42kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a protein that plays a role in the regulation of glomerular permeability. Mutations in this gene cause steroid-resistant nephrotic syndrome. Alternative splicing results in multiple transcript variants.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: PDCN; SRN1; NPHS2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 284-383 of human NPHS2 (Q9NP85).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.