

CAB23024

Product Information

| | | |
|------------------------------|-----------------------------|------------------------------|
| Product SKU: CAB23024 | Gene ID: 2302 | Size: 20uL, 100uL |
| Clone No: - | Host Species: Rabbit | Reactivity: Mouse,Rat |

Additional Information

| | |
|-----------------------------|--------------------------------|
| Observed MW: 55kDa/ | Conjugate: Unconjugated |
| Calculated MW: 45kDa | Isotype: IgG |

Immunogen Information

Background: This gene encodes a member of the forkhead family of transcription factors. Similar genes in zebrafish and mouse have been shown to regulate the transcription of genes that control the production of motile cilia. The mouse ortholog also functions in the determination of left-right asymmetry. Polymorphisms in this gene are associated with systemic lupus erythematosus and allergic rhinitis.

Recommended Dilution: WB,1:1000 - 1:5000 IP,0.5µg-4µg antibody for 400µg-600µg extracts of whole cells

Synonyms: HFH4; HFH-4; CILD43; FKHL13; FOXJ1

Purification Method: Affinity purification

Immunogen: Recombinant protein of human FOXJ1.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.