

CAB23045

Product Information

Product SKU: CAB23045	Gene ID: 58205	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Mouse,Rat

Additional Information

Observed MW: 45-60kDa	Conjugate: Unconjugated
Calculated MW: 28kDa	Isotype: IgG

Immunogen Information

Background: Acts upstream of or within negative regulation of T cell proliferation and positive regulation of T cell proliferation. Predicted to be located in plasma membrane. Predicted to be integral component of membrane. Predicted to be active in external side of plasma membrane. Is expressed in several structures, including alimentary system epithelium; forebrain; nose; retina; and skin. Orthologous to human PDCD1LG2 (programmed cell death 1 ligand 2).

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: Btdc; B7-DC; PD-L2; F730015O22Rik; PD-L2/CD273/B7-DC

Purification Method: Affinity purification

Immunogen: Recombinant protein of mouse PD-L2/CD273/B7-DC.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.