

CAB23047

Product Information

Product SKU: CAB23047	Gene ID: 55074	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 139kDa	Conjugate: Unconjugated
Calculated MW: 98kDa	Isotype: IgG

Immunogen Information

Background: Predicted to enable oxidoreductase activity. Predicted to be involved in response to oxidative stress. Predicted to act upstream of or within several processes, including adult walking behavior; negative regulation of neuron death; and negative regulation of peptidyl-cysteine S-nitrosylation. Predicted to be located in mitochondrion and nucleolus. Predicted to be active in nucleus. Implicated in cerebellar hyplasia/atrophy, epilepsy, and global developmental delay.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: TLDC3; CHEGDD; Nbla00307; OXR1

Purification Method: Affinity purification

Immunogen: Recombinant protein of human OXR1.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.