

CAB23048

Product Information

Product SKU:	CAB23048	Gene ID:	51303	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	22kDa	Conjugate:	Unconjugated
Calculated MW:	98kDa	Isotype:	IgG

Immunogen Information

Background: FKBP11 belongs to the FKBP family of peptidyl-prolyl cis/trans isomerases, which catalyze the folding of proline-containing polypeptides. The peptidyl-prolyl isomerase activity of FKBP proteins is inhibited by the immunosuppressant compounds FK506 and rapamycin (Rulten et al., 2006 [PubMed 16596453]).

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: FKBP19; FKBP11

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 28-150 of human FKBP11 (NP_057678.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.