

CAB23049

Product Information

Product SKU:	CAB23049	Gene ID:	19124	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	45kDa	Conjugate:	Unconjugated
Calculated MW:	27kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable protein C-terminus binding activity. Predicted to be involved in negative regulation of coagulation. Predicted to act upstream of or within blood coagulation. Located in centrosome. Is expressed in several structures, including blastocyst; endocardial tissue; female reproductive system; liver; and vascular element. Orthologous to human PROCR (protein C receptor).

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: Ccca; Epcr; Ccd41; CD201

Purification Method: Affinity purification

Immunogen: Recombinant protein of mouse CD201.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.