

CAB23050

Product Information

Product SKU: CAB23050	Gene ID: 17470	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Mouse,Rat

Additional Information

Observed MW: 45-50kDa	Conjugate: Unconjugated
Calculated MW: 31kDa	Isotype: IgG

Immunogen Information

Background: Predicted to enable protein binding activity involved in heterotypic cell-cell adhesion. Involved in negative regulation of cell population proliferation; negative regulation of macrophage activation; and negative regulation of neuroinflammatory response. Located in cell surface. Is expressed in several structures, including anterior naris; aorta; bone; nervous system; and retina. Used to study autoimmune disease. Orthologous to human CD200 (CD200 molecule).

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: OX2; Mox2; CD200

Purification Method: Affinity purification

Immunogen: Recombinant protein of mouse CD200.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.