

CAB23053

Product Information

Product SKU:	CAB23053	Gene ID:	16407	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	150kDa	Conjugate:	Unconjugated
Calculated MW:	129kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable metal ion binding activity. Predicted to act upstream of or within cell adhesion and integrin-mediated signaling pathway. Located in external side of plasma membrane. Is expressed in several structures, including alimentary system; genitourinary system; lymph node; nasal septum; and nucleus pulposus. Orthologous to human ITGAE (integrin subunit alpha E).

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: CD103; aM290; alpha-E1; A530055J10; alpha-M290

Purification Method: Affinity purification

Immunogen: Recombinant protein of mouse CD103.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.