

CAB23059

Product Information

Product SKU: CAB23059	Gene ID: 2358	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: Refer to figures	Conjugate: Unconjugated
Calculated MW: 39kDa	Isotype: IgG

Immunogen Information

Background: Enables amyloid-beta binding activity; scavenger receptor binding activity; and signaling receptor activity. Involved in several processes, including cellular response to amyloid-beta; positive regulation of monocyte chemotaxis; and regulation of defense response. Located in cytoplasm and plasma membrane.

Recommended Dilution: IF/ICC, 1:50 - 1:200

Synonyms: ALX; ALXR; HM63; FMLPX; FPR2A; FPRH1; FPRH2; FPRL1; LXA4R; FMLP-R-II; FPR2

Purification Method: Affinity purification

Immunogen: Recombinant protein of human FPR2.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.