

## CAB23063

---

### Product Information

<b>Product SKU:</b>	CAB23063	<b>Gene ID:</b>	197259	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	54kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	54kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

<b>Background:</b>	This gene belongs to the protein kinase superfamily. The encoded protein contains a protein kinase-like domain; however, is thought to be inactive because it lacks several residues required for activity. This protein plays a critical role in tumor necrosis factor (TNF)-induced necroptosis, a programmed cell death process, via interaction with receptor-interacting protein 3 (RIP3), which is a key signaling molecule in necroptosis pathway. Inhibitor studies and knockdown of this gene inhibited TNF-induced necrosis. High levels of this protein and RIP3 are associated with inflammatory bowel disease in children. Alternatively spliced transcript variants have been described for this gene.
<b>Recommended Dilution:</b>	WB, 1:500 - 1:1000
<b>Synonyms:</b>	hMLKL; KL
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 179-471 of human MLKL. (NP_689862.1).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.