

CAB23074

Product Information

Product SKU: CAB23074	Gene ID: -	Size: 20uL, 100uL
Clone No: ARC56943	Host Species: Rabbit	Reactivity: Goat

Additional Information

Observed MW: 62kDa	Conjugate: Unconjugated
Calculated MW: -	Isotype: IgG

Immunogen Information

Background:

Recommended Dilution: WB, 1:2000 - 1:6000

Synonyms: -

Purification Method: Affinity purification

Immunogen: Goat IgG

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.