

CAB23078

Product Information

Product SKU:	CAB23078	Gene ID:	8694	Size:	20uL, 100uL
Clone No:	ARC58696	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	55kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a multipass transmembrane protein that functions as a key metabolic enzyme. The encoded protein catalyzes the conversion of diacylglycerol and fatty acyl CoA to triacylglycerol. This enzyme can also transfer acyl CoA to retinol. Activity of this protein may be associated with obesity and other metabolic diseases.

Recommended Dilution: IHC-P,1:50 - 1:200

Synonyms: ARAT; DGAT; ARGPI; DIAR7; DGAT1

Purification Method: Affinity purification

Immunogen: Recombinant protein of human DGAT1.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.