

CAB23087

Product Information

Product SKU: CAB23087	Gene ID: 929017	Size: 20uL, 100uL
Clone No: ARC59811	Host Species: Rabbit	Reactivity: Monkeypox virus

Additional Information

Observed MW: 40-50kDa	Conjugate: Unconjugated
Calculated MW: 35kDa	Isotype: IgG

Immunogen Information

Background:

Recommended Dilution: WB, 1:2000 - 1:6000

Synonyms: -

Purification Method: Affinity purification

Immunogen: Recombinant protein of human Monkeypox virus (MPXV) E8L.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH 7.3.