

**CAB23089**

---

## Product Information

<b>Product SKU:</b>	CAB23089	<b>Gene ID:</b>	23208	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC59135	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	48kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	48kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** This gene is a member of the synaptotagmin gene family and encodes a protein similar to other family members that are known calcium sensors and mediate calcium-dependent regulation of membrane trafficking in synaptic transmission. The encoded protein is also a substrate for ubiquitin-E3-ligase parkin. The gene has previously been referred to as synaptotagmin XII but has been renamed synaptotagmin XI to be consistent with mouse and rat official nomenclature.

**Recommended Dilution:** WB,1:1000 - 1:5000

**Synonyms:** SYT12; sytXI; Synaptotagmin XI (SYT11)

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of human Synaptotagmin XI (SYT11).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.