

CAB23126

Product Information

| | | | | | |
|---------------------|----------|----------------------|--------|--------------------|-------------|
| Product SKU: | CAB23126 | Gene ID: | 2118 | Size: | 20uL, 100uL |
| Clone No: | ARC58618 | Host Species: | Rabbit | Reactivity: | Human |

Additional Information

| | | | |
|-----------------------|-------|-------------------|--------------|
| Observed MW: | 61kDa | Conjugate: | Unconjugated |
| Calculated MW: | 54kDa | Isotype: | IgG |

Immunogen Information

Background: Enables DNA-binding transcription activator activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Involved in positive regulation of keratinocyte differentiation and positive regulation of transcription by RNA polymerase II. Located in chromosome and nucleolus.

Recommended Dilution: WB, 1:2000 - 1:9000

Synonyms: E1AF; PEA3; E1A-F; PEAS3; ETV4

Purification Method: Affinity purification

Immunogen: A synthesized peptide derived from human ETV4.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.