

## CAB23148

---

### Product Information

<b>Product SKU:</b>	CAB23148	<b>Gene ID:</b>	22925	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC59749	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	182kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	169kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene represents a phospholipase A2 receptor. The encoded protein likely exists as both a transmembrane form and a soluble form. The transmembrane receptor may play a role in clearance of phospholipase A2, thereby inhibiting its action. Polymorphisms at this locus have been associated with susceptibility to idiopathic membranous nephropathy. Alternatively spliced transcript variants encoding different isoforms have been identified.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** PLA2R; PLA2-R; PLA2IR; CLEC13C; PLA2G1R; PLA2R1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of human PLA2R1.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.