

## CAB23162

---

### Product Information

<b>Product SKU:</b>	CAB23162	<b>Gene ID:</b>	80312	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC59030	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

---

### Additional Information

<b>Observed MW:</b>	310kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	235kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** DNA methylation is an epigenetic mechanism that is important for controlling gene expression. The protein encoded by this gene is a demethylase that belongs to the TET (ten-eleven translocation) family. Members of the TET protein family play a role in the DNA methylation process and gene activation.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** LCX; CXXC6; bA119F7.1; TET1

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 1800-1900 of human TET (NP\_085128.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.