

CAB23163

Product Information

Product SKU:	CAB23163	Gene ID:	7424	Size:	20uL, 100uL
Clone No:	ARC59857	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	47kDa	Conjugate:	Unconjugated
Calculated MW:	47kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a member of the platelet-derived growth factor/vascular endothelial growth factor (PDGF/VEGF) family. The encoded protein promotes angiogenesis and endothelial cell growth, and can also affect the permeability of blood vessels. The proprotein is further cleaved into a fully processed form that can bind and activate VEGFR-2 and VEGFR-3 receptors.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: VRP; Flt4-L; LMPH1D; LMPHM4; VEGFC

Purification Method: Affinity purification

Immunogen: A synthesized peptide derived from human VEGFC.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.