

CAB23240

Product Information

Product SKU:	CAB23240	Gene ID:	1105	Size:	20uL, 100uL
Clone No:	ARC59018	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	250kDa	Conjugate:	Unconjugated
Calculated MW:	197kDa	Isotype:	IgG

Immunogen Information

Background: The CHD family of proteins is characterized by the presence of chromo (chromatin organization modifier) domains and SNF2-related helicase/ATPase domains. CHD genes alter gene expression possibly by modification of chromatin structure thus altering access of the transcriptional apparatus to its chromosomal DNA template.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: CHD-1; PILBOS; CHD1

Purification Method: Affinity purification

Immunogen: Recombinant protein of human CHD1.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.