

CAB23252

Product Information

| | | | | | |
|---------------------|----------|----------------------|--------|--------------------|-------------|
| Product SKU: | CAB23252 | Gene ID: | 797 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Mouse |

Additional Information

| | | | |
|-----------------------|-------|-------------------|--------------|
| Observed MW: | 16kDa | Conjugate: | Unconjugated |
| Calculated MW: | 13kDa | Isotype: | IgG |

Immunogen Information

Background: Predicted to enable calcitonin receptor binding activity. Predicted to be involved in adenylate cyclase-activating G protein-coupled receptor signaling pathway and regulation of cytosolic calcium ion concentration. Predicted to be located in extracellular region. Predicted to be active in extracellular space.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: CALC2; CGRP2; CGRP-II

Purification Method: Affinity purification

Immunogen: Recombinant protein of Human CGPR.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.