

CAB23253

Product Information

Product SKU:	CAB23253	Gene ID:	16783	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	90-120kDa	Conjugate:	Unconjugated
Calculated MW:	44kDa	Isotype:	IgG

Immunogen Information

Background: Enables protein domain specific binding activity. Involved in protein stabilization. Located in several cellular components, including cytoplasmic vesicle; sarcolemma; and vacuole. Is expressed in several structures, including 1-cell stage embryo; central nervous system; craniocervical region bone; extraembryonic component; and gut. Orthologous to human LAMP1 (lysosomal associated membrane protein 1).

Recommended Dilution: WB,1:100 - 1:500 IF/ICC,1:50 - 1:200

Synonyms: P2B; Perk; LGP-A; CD107a; Lamp-1; LGP-120; LAMP1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 208-363 of mouse LAMP1 (NP_001304282.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.