

CAB23269

Product Information

Product SKU:	CAB23269	Gene ID:	284207	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	26kDa/34kDa	Conjugate:	Unconjugated
Calculated MW:	34kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable hormone activity. Predicted to be involved in several processes, including brown fat cell differentiation; energy homeostasis; and positive regulation of brown fat cell differentiation. Located in extracellular exosome.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: -

Purification Method: Affinity purification

Immunogen: Recombinant protein of mouse METRNL.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.