

CAB23276

Product Information

Product SKU:	CAB23276	Gene ID:	8505	Size:	20uL, 100uL
Clone No:	ARC59258	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	111kDa	Isotype:	IgG

Immunogen Information

Background: Poly(ADP-ribose) glycohydrolase (PARG) is the major enzyme responsible for the catabolism of poly(ADP-ribose), a reversible covalent-modifier of chromosomal proteins. The protein is found in many tissues and may be subject to proteolysis generating smaller, active products. Several transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: IHC-P,1:50 - 1:200

Synonyms: PARG99; PARG?

Purification Method: Affinity purification

Immunogen: Recombinant protein of human PARG.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.