

## CAB23299

---

### Product Information

<b>Product SKU:</b>	CAB23299	<b>Gene ID:</b>	16590	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC60558	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	120kDa/145kDa	<b>Conjugate:</b>	-
<b>Calculated MW:</b>	23KDa/109KDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** The c-Kit proto-oncogene is the cellular homolog of the transforming gene of a feline retrovirus (v-Kit). The c-kit protein includes characteristics of a protein kinase transmembrane receptor. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** W; Bs; Fdc; Ssm; SCO1; SCO5; SOW3; CD117; c-KIT; Tr-kit; Gsfsc01; Gsfsc05; Gsfscow3; CD117/c-Kit

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of mouse CD117/c-Kit.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.