

CAB23428

Product Information

Product SKU:	CAB23428	Gene ID:	3693	Size:	20uL, 100uL
Clone No:	ARC60869	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	100kDa	Conjugate:	Unconjugated
Calculated MW:	88kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a beta subunit of integrin, which can combine with different alpha chains to form a variety of integrin heterodimers. Integrins are integral cell-surface receptors that participate in cell adhesion as well as cell-surface mediated signaling. The α v β 5 integrin is involved in adhesion to vitronectin.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: ITGB5; integrin beta-5; Integrin β 5/ITGB5

Purification Method: Affinity purification

Immunogen: Recombinant protein of human Integrin β 5/ITGB5.

Storage: Store at -20°C . Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH 7.3.