

CAB23447

Product Information

Product SKU:	CAB23447	Gene ID:	6656	Size:	20uL, 100uL
Clone No:	ARC57513	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	39kDa	Isotype:	IgG

Immunogen Information

Background: This intronless gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional activator after forming a protein complex with other proteins. In mice, a similar protein regulates the gamma-crystallin genes and is essential for lens development.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: SOX1; SRY-box 1

Purification Method: Affinity purification

Immunogen: Recombinant protein of human SOX1.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.