

CAB23499

Product Information

Product SKU:	CAB23499	Gene ID:	2549	Size:	20uL, 100uL
Clone No:	ARC52035	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	110kDa	Conjugate:	Unconjugated
Calculated MW:	76kDa/80kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a member of the IRS1-like multisubstrate docking protein family. It is an important mediator of branching tubulogenesis and plays a central role in cellular growth response, transformation and apoptosis. Two transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: DFNB26; GAB1

Purification Method: Affinity purification

Immunogen: Recombinant protein of human GAB1.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.