

CAB23539

Product Information

Product SKU:	CAB23539	Gene ID:	64863	Size:	20uL, 100uL
Clone No:	ARC59459	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	54kDa	Isotype:	IgG

Immunogen Information

Background: Enables RNA methyltransferase activity and site-specific DNA-methyltransferase (adenine-specific) activity. Involved in nucleic acid metabolic process; regulation of RNA metabolic process; and regulation of mitochondrial DNA replication. Located in cytosol; mitochondrial matrix; and nucleus.

Recommended Dilution: IHC-P,1:50 - 1:200

Synonyms: HsT661; METTL4

Purification Method: Affinity purification

Immunogen: Recombinant protein of human METTL4.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.