

CAB23542

Product Information

Product SKU:	CAB23542	Gene ID:	-	Size:	20uL, 100uL
Clone No:	ARC60394	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	60kDa	Conjugate:	Unconjugated
Calculated MW:	37KDa/42KDa	Isotype:	IgG

Immunogen Information

Background:

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: -

Purification Method: Affinity purification

Immunogen: Recombinant protein of human Human IgA1.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.