

CAB23553

Product Information

Product SKU:	CAB23553	Gene ID:	51155	Size:	20uL, 100uL
Clone No:	ARC3106	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	25kDa	Conjugate:	Unconjugated
Calculated MW:	16kDa	Isotype:	IgG

Immunogen Information

Background: Located in nuclear membrane; nucleolus; and nucleoplasm.

Recommended Dilution: WB,1:200 - 1:500

Synonyms: HN1; ARM2; HN1A

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 50-154 of human HN1 (NP_057269.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.