

CAB23655

Product Information

Product SKU:	CAB23655	Gene ID:	17086	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	46kDa	Conjugate:	Unconjugated
Calculated MW:	37kDa	Isotype:	IgG

Immunogen Information

Background: Acts upstream of or within defense response to virus and detection of virus. Located in cell surface. Orthologous to human NCR1 (natural cytotoxicity triggering receptor 1).

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: Ly94; NKp46; CD335

Purification Method: Affinity purification

Immunogen: Recombinant protein of mouse CD335.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.