

CAB23660

Product Information

Product SKU:	CAB23660	Gene ID:	25375	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	55kDa	Conjugate:	Unconjugated
Calculated MW:	50kDa	Isotype:	IgG

Immunogen Information

Background: Enables 3-chloroallyl aldehyde dehydrogenase activity. Involved in several processes, including response to cAMP; response to glucocorticoid; and response to hypoxia. Located in cytosol. Human ortholog(s) of this gene implicated in arteriosclerosis. Orthologous to human ALDH3A1 (aldehyde dehydrogenase 3 family member A1).

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: AHDC; Aldh; Ahd-c; Aldh3; Aldh3a1

Purification Method: Affinity purification

Immunogen: Recombinant protein of mouse Aldh3a1.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.