

## CAB23665

---

### Product Information

|                     |          |                      |        |                    |             |
|---------------------|----------|----------------------|--------|--------------------|-------------|
| <b>Product SKU:</b> | CAB23665 | <b>Gene ID:</b>      | 50929  | <b>Size:</b>       | 20uL, 100uL |
| <b>Clone No:</b>    | -        | <b>Host Species:</b> | Rabbit | <b>Reactivity:</b> | Human,Mouse |

---

### Additional Information

|                       |             |                   |              |
|-----------------------|-------------|-------------------|--------------|
| <b>Observed MW:</b>   | 16kDa/18kDa | <b>Conjugate:</b> | Unconjugated |
| <b>Calculated MW:</b> | 20kDa       | <b>Isotype:</b>   | IgG          |

---

### Immunogen Information

**Background:** Predicted to enable cytokine activity. Acts upstream of or within several processes, including negative regulation of inflammatory response; reactive oxygen species metabolic process; and regulation of tyrosine phosphorylation of STAT protein. Predicted to be located in extracellular space. Is expressed in retina and skin. Used to study liposarcoma. Human ortholog(s) of this gene implicated in asthma. Orthologous to human IL22 (interleukin 22).

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** IL-22; If2b1; Iltif; IL-22a; ILTIFa; IL22

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of mouse IL22.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.