

## CAB23693

---

### Product Information

<b>Product SKU:</b>	CAB23693	<b>Gene ID:</b>	7091	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC61324	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	84kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	84kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Predicted to enable transcription corepressor activity. Predicted to be involved in negative regulation of canonical Wnt signaling pathway. Predicted to act upstream of or within Wnt signaling pathway; cellular response to leukemia inhibitory factor; and negative regulation of transcription by RNA polymerase II. Located in nucleoplasm. Part of beta-catenin-TCF complex.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** ESG; BCE1; ESG4; GRG4; BCE-1; Grg-4; E(spl); E(spl); TLE4

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-138 of human TLE4 (NP\_001269677).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.