

CAB23761

Product Information

Product SKU:	CAB23761	Gene ID:	2804	Size:	20uL, 100uL
Clone No:	ARC61107	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	376KDa	Isotype:	IgG

Immunogen Information

Background: Enables RNA binding activity. Involved in protein localization to pericentriolar material. Located in Golgi apparatus and endoplasmic reticulum-Golgi intermediate compartment.

Recommended Dilution: IHC-P,1:100 - 1:500

Synonyms: GCP; GCP372; GOLIM1; Giantin

Purification Method: Affinity purification

Immunogen: Recombinant protein of human Giantin.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.