

CAB23779

Product Information

Product SKU:	CAB23779	Gene ID:	9398	Size:	20uL, 100uL
Clone No:	ARC62044	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	140kDa	Conjugate:	Unconjugated
Calculated MW:	115kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable hydrolase activity, acting on carbon-nitrogen (but not peptide) bonds, in cyclic amides. Predicted to be involved in cell surface receptor signaling pathway. Predicted to act upstream of or within positive regulation of myeloid leukocyte differentiation. Located in extracellular exosome.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: V7; IGSF2; EWI-101; CD101

Purification Method: Affinity purification

Immunogen: Recombinant protein of human CD101.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.