

CAB23786

Product Information

Product SKU: CAB23786	Gene ID: 4837	Size: 20uL, 100uL
Clone No: ARC61849	Host Species: Rabbit	Reactivity: Human,Rat

Additional Information

Observed MW: 28kDa	Conjugate: Unconjugated
Calculated MW: 29kDa	Isotype: IgG

Immunogen Information

Background: N-methylation is one method by which drug and other xenobiotic compounds are metabolized by the liver. This gene encodes the protein responsible for this enzymatic activity which uses S-adenosyl methionine as the methyl donor.

Recommended Dilution: WB,1:1000 - 1:5000

Synonyms: NNMT; nicotinamide N-methyltransferase

Purification Method: Affinity purification

Immunogen: Recombinant protein of human NNMT.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.