

CAB23792

Product Information

Product SKU:	CAB23792	Gene ID:	51377	Size:	20uL, 100uL
Clone No:	ARC60845	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	37kDa	Conjugate:	Unconjugated
Calculated MW:	36kDa/37kDa	Isotype:	IgG

Immunogen Information

Background: Enables endopeptidase inhibitor activity; proteasome binding activity; and thiol-dependent deubiquitinase. Involved in negative regulation of proteasomal ubiquitin-dependent protein catabolic process; positive regulation of smoothed signaling pathway; and protein deubiquitination. Located in cytosol; nucleolus; and nucleoplasm. Colocalizes with Ino80 complex.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: UCHL5; CGI-70; INO80R; UCH-L5; UCH37; ubiquitin C-terminal hydrolase L5; UCH37/UCHL5

Purification Method: Affinity purification

Immunogen: Recombinant protein of human UCH37/UCHL5.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.