

CAB23885

Product Information

Product SKU:	CAB23885	Gene ID:	17067	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	15kDa	Conjugate:	-
Calculated MW:	14KDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable acetylcholine receptor binding activity and acetylcholine receptor inhibitor activity. Predicted to be involved in acetylcholine receptor signaling pathway. Located in external side of plasma membrane. Is expressed in embryo mesenchyme; gut; liver; metanephros; and nasal cavity. Orthologous to human LY6H (lymphocyte antigen 6 family member H).

Recommended Dilution: WB, 1:1000 - 1:5000

Synonyms: Ly6c; Ly-6C; Ly-6C1

Purification Method: Affinity purification

Immunogen: Recombinant protein of mouse Ly-6C.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.