

CAB2389

Product Information

Product SKU:	CAB2389	Gene ID:	10713	Size:	20uL, 100uL
Clone No:	ARC2588	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	65kDa	Conjugate:	Unmodified
Calculated MW:	65kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable thiol-dependent deubiquitinase and zinc ion binding activity. Involved in spliceosomal complex assembly. Located in nucleoplasm. Part of U4/U6 x U5 tri-snRNP complex.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: 65K; SAD1; CGI-21; HSPC332; SNRNP65; USP39

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 466-565 of Human USP39 (Q53GS9).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.