

CAB2390

Product Information

Product SKU:	CAB2390	Gene ID:	1400	Size:	20uL, 100uL
Clone No:	ARC1918	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	62kDa	Conjugate:	Unconjugated
Calculated MW:	62kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of a family of cytosolic phosphoproteins expressed exclusively in the nervous system. The encoded protein is thought to be a part of the semaphorin signal transduction pathway implicated in semaphorin-induced growth cone collapse during neural development. Alternative splicing results in multiple transcript variants.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: DRP1; DRP-1; CRMP-1; DPYSL1; ULIP-3; CRMP1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 450-550 of human CRMP1 (Q14194).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.