

CAB23911

Product Information

Product SKU:	CAB23911	Gene ID:	92140	Size:	20uL, 100uL
Clone No:	ARC62211	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	64kDa	Isotype:	IgG

Immunogen Information

Background: Enables NF-kappaB binding activity; double-stranded RNA binding activity; and transcription coactivator activity. Involved in several processes, including lipopolysaccharide-mediated signaling pathway; positive regulation of intracellular signal transduction; and regulation of transcription, DNA-templated. Located in endoplasmic reticulum; nuclear lumen; and perinuclear region of cytoplasm. Implicated in hepatocellular carcinoma.

Recommended Dilution: IHC-P,1:50 - 1:200

Synonyms: 3D3; AEG1; AEG-1; LYRIC; LYRIC/3D3; LYRIC/AEG1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 308-407 of human LYRIC/AEG1 (NP_848927.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.