

CAB23956

Product Information

Product SKU:	CAB23956	Gene ID:	16181	Size:	20uL, 100uL
Clone No:	ARC62499	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	19kDa	Conjugate:	Unconjugated
Calculated MW:	20kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable cytokine receptor binding activity. Acts upstream of or within insulin secretion and lipid metabolic process. Located in extracellular region and vesicle. Human ortholog(s) of this gene implicated in several diseases, including autoimmune disease (multiple); eye disease (multiple); kidney failure (multiple); lung disease (multiple); and vascular disease (multiple). Orthologous to human IL1RN (interleukin 1 receptor antagonist).

Recommended Dilution: WB,1:1000 - 1:5000

Synonyms: IL-1ra; F630041P17Rik; IL-1RA

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 27-178 of mouse IL-1RA (NP_001034790.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.