

CAB23959

Product Information

Product SKU:	CAB23959	Gene ID:	-	Size:	20uL, 100uL
Clone No:	ARC62861	Host Species:	Rabbit	Reactivity:	Widely Range,Oryza sativa,Arabidopsis thaliana

Additional Information

Observed MW:	42kDa	Conjugate:	Unconjugated
Calculated MW:	42kDa	Isotype:	IgG

Immunogen Information

Background: Actin is an ancient conserved protein that is encoded by multiple isoforms in multi-cellular organisms. There are eight functional actin genes in the Arabidopsis genome.

Recommended Dilution: WB,1:2000 - 1:10000

Synonyms: ACT12; Actin(Plant Specific)

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of Actin(Plant Specific)(NP_190236.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.