

**CAB23994**

---

**Product Information**

<b>Product SKU:</b>	CAB23994	<b>Gene ID:</b>	51131	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC62550	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

**Additional Information**

<b>Observed MW:</b>	34kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	38kDa	<b>Isotype:</b>	IgG

---

**Immunogen Information**

**Background:** This gene encodes a protein containing a PHD (plant homeodomain) type zinc finger. This gene has been identified in some studies as a candidate gene for asthma. Naturally-occurring readthrough transcription may occur from the upstream SETDB2 (SET domain bifurcated 2) gene to this locus. Alternative splicing results in multiple transcript variants.

**Recommended Dilution:** WB,1:1000 - 1:5000

**Synonyms:** APY; BCAP; IGEL; IGER; IGHF; NYREN34; NY-REN-34; [KD Validated] PHF11

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 40-321 of human PHF11 (NP\_001035533.1) .

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.