

CAB24021

Product Information

Product SKU:	CAB24021	Gene ID:	10241	Size:	20uL, 100uL
Clone No:	ARC54603	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	60kDa	Conjugate:	Unconjugated
Calculated MW:	52kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a coiled-coil domain-containing protein. The encoded protein functions as a receptor for ubiquitin-coated bacteria and plays an important role in innate immunity by mediating macroautophagy. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

Recommended Dilution: WB, 1:1000 - 1:5000

Synonyms: NDP52

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human NDP52 (NP_005822.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH 7.3.