

CAB24043

Product Information

Product SKU: CAB24043	Gene ID: 23428	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 72kDa	Conjugate: Unconjugated
Calculated MW: 58kDa	Isotype: IgG

Immunogen Information

Background: Enables several functions, including neutral amino acid transmembrane transporter activity; thyroid hormone transmembrane transporter activity; and toxin transmembrane transporter activity. Involved in L-alanine import across plasma membrane; L-leucine import across plasma membrane; and thyroid hormone transport. Located in plasma membrane. Part of basolateral plasma membrane and microvillus membrane.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: LAT2; LPI-PC1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 201-300 of human SLC7A8 (NP_036376.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.