

CAB24059

Product Information

Product SKU:	CAB24059	Gene ID:	79792	Size:	20uL, 100uL
Clone No:	ARC57989	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	31kDa	Conjugate:	Unconjugated
Calculated MW:	53kDa	Isotype:	IgG

Immunogen Information

Background: Gasdermin D is a member of the gasdermin family. Members of this family appear to play a role in regulation of epithelial proliferation. Gasdermin D has been suggested to act as a tumor suppressor. Alternatively spliced transcript variants have been described.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: DF5L; DFNA5L; FKSG10; GSDMDC1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to residues surrounding Asp275 of human Gasdermin D(NP_079012.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.