

CAB24353

Product Information

Product SKU:	CAB24353	Gene ID:	2131	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	86kDa	Conjugate:	Unconjugated
Calculated MW:	86kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes an endoplasmic reticulum-resident type II transmembrane glycosyltransferase involved in the chain elongation step of heparan sulfate biosynthesis. Mutations in this gene cause the type I form of multiple exostoses.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: EXT; LGS; TTV; LGCR; TRPS2; EXT1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 300-400 of human EXT1(NP_000118.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.