

## CAB24437

---

### Product Information

<b>Product SKU:</b> CAB24437	<b>Gene ID:</b> 4312	<b>Size:</b> 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human,Mouse

---

### Additional Information

<b>Observed MW:</b> 57kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 54kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** This gene encodes a member of the peptidase M10 family of matrix metalloproteinases (MMPs). Proteins in this family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. The encoded preproprotein is proteolytically processed to generate the mature protease. This secreted protease breaks down the interstitial collagens, including types I, II, and III. The gene is part of a cluster of MMP genes on chromosome 11. Mutations in this gene are associated with chronic obstructive pulmonary disease (COPD). Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that is proteolytically processed.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** CLG; CLGN; MMP1

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 369-469 of human MMP1 (NP\_002412.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.